


Patient clinical status over 31 days from Sep 17 to Oct 17

■ HR ■ Breathing



Last 24 Hours

Last 24 Hours: NORMAL

30-Day Tier: YELLOW

| Vital      | Avg     | Min     | Max     |
|------------|---------|---------|---------|
| Heart Rate | 57 bpm  | 50 bpm  | 85 bpm  |
| Breathing  | 16 brpm | 10 brpm | 46 brpm |

Last 24 hours: one elevated breathing episode in the morning, otherwise within normal range.

Episodic Events (last 24h)

| Date   | Time Span     | Duration | Episodes | Condition          | Avg | Min | Max | Comment                                                                  |
|--------|---------------|----------|----------|--------------------|-----|-----|-----|--------------------------------------------------------------------------|
| Oct 17 | 9 AM to 10 AM | 1h       | 1        | Elevated Breathing | 28  | 13  | 46  | Check SpO2 and lung sounds; assess for fluid overload or early infection |

43 episodic events Detected, spanning 50h total. Conditions: Elevated Breathing, Elevated Heart Rate, High Breathing, High Heart Rate, Very High Breathing. **Priority grade: Medium**

| Episodic Events | Total Hours | Highest Burden | Episodes/day | Coverage |
|-----------------|-------------|----------------|--------------|----------|
| 43              | 50h         | High Breathing | 1            | 100%     |

Major Findings: Sustained high breathing (avg 35, peak 58)

| Date   | Time Span    | Duration | Episodes | Condition                   | Avg | Min | Max | Comment                                                                  |
|--------|--------------|----------|----------|-----------------------------|-----|-----|-----|--------------------------------------------------------------------------|
| Sep 24 | 1 PM to 2 PM | 15h      | 12       | High Breathing              | 35  | 12  | 58  | Check SpO2, work of breathing, fever; consider pulmonary cause or CHF    |
| Oct 05 | 2 AM to 3 AM | 3h       | 2        | High Heart Rate             | 104 | 87  | 115 | Check temp, hydration, pain; review for arrhythmia or infection          |
| Oct 08 | 8 AM to 9 AM | 3h       | 3        | High Breathing              | 35  | 12  | 54  | Check SpO2, work of breathing, fever; consider pulmonary cause or CHF    |
| Sep 25 | 7 AM to 8 AM | 6h       | 6        | Elevated Breathing          | 26  | 14  | 57  | Check SpO2 and lung sounds; assess for fluid overload or early infection |
| Oct 05 | 1 AM to 2 AM | 2h       | 2        | Elevated Heart Rate         | 98  | 57  | 108 | Check temp, hydration, pain; review for arrhythmia or infection          |
| Oct 06 | 1 PM to 2 PM | 9h       | 8        | Elevated Breathing          | 26  | 13  | 46  | Check SpO2 and lung sounds; assess for fluid overload or early infection |
| Oct 08 | 7 AM to 8 AM | 2h       | 2        | Very High Breathing (brief) | 46  | 14  | 54  | Check SpO2, work of breathing, fever; consider pulmonary cause or CHF    |


The P5 to P95 heart rate spread was 33 bpm (52 to 85), indicating significant cardiac variability.

Suggested Clinical Review Actions

- High Breathing (avg 35 brpm, peak 58 brpm): Check SpO2, work of breathing, fever; consider pulmonary cause or CHF
- Elevated Breathing (avg 26 brpm, peak 46 brpm): Check SpO2 and lung sounds; assess for fluid overload or early infection
- High Heart Rate (avg 104 bpm, peak 115 bpm): Check temp, hydration, pain; review for arrhythmia or infection
- Elevated Heart Rate (avg 98 bpm, peak 108 bpm): Check temp, hydration, pain; review for arrhythmia or infection
- Very High Breathing (avg 46 brpm, peak 54 brpm): Check SpO2, work of breathing, fever; consider pulmonary cause or CHF


Clinical Pattern Observations:

- 1. **High variability:** HR spread 32 bpm (52 to 85).
- 2. **Clustered pattern:** Episodic events on 5 of 31 days (16%).
- 3. **Daytime pattern:** Episodes concentrated during waking hours.




Red bars indicate days with detected episodic events.

Vital Sign Distribution



Monitoring Activity



Bars indicate monitoring hours/day. Red threshold line indicates clinical monitoring baseline.

■ HR episode ■ RR episode ■ Recorded, no episode ■ No data

Clinical Alerting Thresholds

- Very Low HR < 40 bpm
- Low HR < 45 bpm
- Elevated HR > 95 bpm
- High HR > 100 bpm
- Very High HR > 110 bpm
- Elevated Breathing > 24 brpm
- High Breathing > 30 brpm
- Very High Breathing > 40 brpm

Episodic event: threshold exceeded continuously for  $\geq 1$  hour. Episodes separated by  $\leq 1$  hour of normal vitals are counted as the same event. Brief threshold crossings that do not sustain are not flagged.